

## Form to release a Biosafety Cabinet (BSC)

| GENERAL INFORMATION   |  |
|---|--|
| Date and time of repair   |  |
| Department  |  |
| Room number   |  |
| Containment level *   |  |
| BSC registration number   |  |
| Description fault/maintenance   |  |
| DECLARATION   |  |
| <p>The parts of the BSC that may be contaminated (with GMOs)**, have been disinfected according to the requirements specified in the work instruction <a href="#">Disinfection of laboratories</a>.</p> <p>The disinfection has been carried out with the following disinfectant (concentration/name disinfectant):</p> <p>.....</p> <p>Therefore, no microbiological risks are associated with repairing the BSC. Please consult the Biological Safety Officer if you have any doubts about the level of disinfection when opening the BSC.</p> <p><b>In addition, I declare that:</b></p> <p><input type="checkbox"/> the pre-filter is checked and <b>doesn't need</b> to be replaced or</p> <p><input type="checkbox"/> the pre-filter is checked and <b>needs</b> to be replaced</p> |  |
| SIGNATURE ABV   |  |
| Name  |  |
| Department  |  |
| Date/time of release  |  |
| Signature   |  |
| SIGNATURE BSO   |  |
| Name  |  |
| Date/time of release  |  |
| Signature   |  |

\* For containment level I, the ABV signs; for containment level II, the BSO and the ABV sign

\*\* Top and bottom of work surface, (integrated) drip tray around the pre-filter